

1. Part I: Relaxation time

Pazis Huth, Zinich Figueroa

In this part of the experiment we will measure the relaxation time T_1 & T_2 for two samples: GDSO¹, GDSO². For that we use a different software that (automat) to perform 10 measurements. and The mean and standard deviation will be stored in a txt file.

~~The measurements are then for~~

For the measurement of T_2 we will use two methods: spin-echo and Carr-Purcell sequence. In both case the time between the 90° & 180° τ to sample the free induction decay curve.

At every measurement we keep the working frequency $\nu_w \approx 1 \text{ kHz} \pm 50 \text{ kHz}$

The final result are written in Tab 1.

Tab 1: relaxation time T_1 & T_2 for GDS samples

	$T_{2sp} [\text{ms}]$	$T_{2CP} [\text{ms}]$	$T_1 [\text{ms}]$	
GDSO ¹	104 ± 1	$102,8 \pm 2,1$	$159,9 \pm 1,3$	$T_{2CO} = (111,3 \pm 1,5) \text{ ms}$
GDSO ²	$115,8 \pm 1,2$	$116,9 \pm 0,9$	$154,4 \pm 1,2$	$T_{2SP} = (110,1 \pm 6) \text{ ms}$

probably due to some the temperature susceptibility the measurement was wrong.

2. Part II: Chemical shift

In this part we want to use the characteristic chemical shift of substances to identify ~~#~~ ^S samples. For every sample we also have a second one with a reference substance, in our case TMS.

To identify the substances we measure the spin echo and perform a Fourier transform. In order to ensure a clear measurement we will use a air pump to make the sample rotate, in this way we minimize the inhomogeneities of the magnetic field.

(Then) We then fit a gauss and measure the position

of the peak. ~~The~~ Using the By comparing the measurement of the sample with the reference and without the reference we can identify the peak corresponding to the reference substance. we then measure the difference from the ^{other} peaks of and reference. Using this data we then use the Fig. 12 of the script to identify the active group:

substance C: $\Delta \text{ppm}_1 = 2.2$; $\Delta \text{ppm}_2 = 11.8$ } C $\hat{=}$ acetic acid
~~at~~ $\hookrightarrow \text{CH}_3$ $\hookrightarrow \text{COOH}$

substance D

we measure two peaks $\rightarrow \Delta \text{ppm}_1 = 4.1$

since) and \Rightarrow ~~the~~ fluoroacetonitril $\Delta \text{ppm}_2 = 6.5$

substance E: $\Delta \text{ppm}_1 = 2.3$; $\Delta \text{ppm}_2 = 7.1$

substance A:

we see three peaks at the fourier transform without reference, therefore we can ~~immediately~~ immediately recognize it as fluoroacetone.

substance B:

$\Delta \text{ppm}_1 = 2.5$; $\Delta \text{ppm}_2 = 7.2$

\Rightarrow ~~to~~ we find the same spectrum for substance ~~A~~ B, however we observe a higher ~~an~~ intensity for the ~~reverse~~ first peak, therefore ~~the~~ in the spectrum of B therefore ~~the~~ we identify it with p-xylol and E $\hat{=}$ toluol.

Part III: Imaging

1. In this experiment we will start by taking a 1-dim profile of ~~the~~ three samples. ~~At~~ Two of them contains the same oil, but with a higher volume. We use this two to check the linearity range of our machine. *

With these comparison we ~~for~~ estimate the linearity range. $[15 \text{ mm}, 15 \text{ m}]$.

Then we insert a third sample with which is a teflon immersed in oil. Since teflon does not emit a signal we expect a profile made up of valley where the teflon is and peaks where the oil is at.

* * as expected we observe a (constant) a square signal with some noise, but for fairly constant.

In the other hand the large sample ~~to~~ ~~it~~ show an increase in the intensity of the external signals, which shows the deviation from linearity.

(per collision)

diffusion

$$\frac{\partial C}{\partial t} = \Delta \cdot \frac{\partial^2 C}{\partial x^2}$$

\sim
 > 0

2. Secondly we prepare a glass tube with diam 1.5 cm of sand and the pour oil in it. ~~we~~ we then proceed to insert it in our imaging device to ~~it~~ study the time evolution.

of the position

After a while we can observe that the curve ^{of the position} is concave, which means that the system's ~~evolution~~ evolution is not a diffusion process.

3.

In the next part of the experiment we will take 2-d profiles of different object.