

F61: Nuclear Magnetic Resonance

Paris J. Huth & Q'inich Figueroa

May 2024

Abstract

In Protocoll we will examine the usage nuclear magnetic resonance to identify probes and reveal the structure of objects.

1 Basics

1.1 Basics of Nuclear Magnetic Resonance

Nuclear Magnetic Resonance techniques relay of the interaction between the magnetic dipole moment

$$\vec{\mu} = \hbar\gamma\vec{S} \quad (1)$$

of nuclei with non-zero spin S and an external magnetic field \vec{B}_0 . In the following paper γ represents the gyromagnetic ratio of protons:

$$\gamma_{proton} = 2.6752 \cdot 10^8 \text{ sec}^{-1} \text{ Tesla}^{-1}. \quad (2)$$

The resulting interaction energy is defines as:

$$\Delta E = -\vec{\mu} \cdot \vec{B}_0. \quad (3)$$

In a classical description, this interaction yields two states for the orientation of the protons's magnetic dipole in the external magnetic field: μ_+ (parallel) and μ_- (antiparalle). For a macroscopic sample of N protons, both numbers of occupied states N_+ and N_- , the sum of which comprises N , can be approximated by a Boltzmann distribution:

$$N_{\pm} = N_0 e^{-\frac{E_0 \pm \Delta E}{kt}} \quad (4)$$

with N_0 as a normalization factor. However $N_+ > N_-$, since the parallel state is energetically favorable. The predominance of protons in the parallel

state leads to a macroscopic magnetization, whose ground state is

$$\vec{M}_0 = \frac{\mu N}{V} \sinh\left(\frac{\mu B}{kT}\right) \vec{e}_z. \quad (5)$$

In our case, a weak field ($\mu B \gg kT$), the former expression simplifies to

$$\vec{M}_0 = \frac{N}{V} \frac{\hbar^2 \gamma^2 I(I+1)}{3kT} \vec{B}_0 \propto \frac{\vec{B}_0}{T}, \quad (6)$$

i.e the law of Curie.

In general, the magnetization can have an arbitrary orientation related to the external magnetic field, however such a system will decay asymptotically into the ground state, since \vec{M}_0 minimizes the energy.

The interaction between the macroscopic magnetization and the external magnetic field result in a torque

$$\vec{\tau} = \gamma \vec{M} \wedge \vec{B}_0. \quad (7)$$

If the magnetization is separated into parallel \vec{M}_\parallel and perpendicular \vec{M}_\perp components, relative to the external field, we quickly see that only the later gives a none trivial expression. Without any relaxation processes, the rate of change of the former is given by

$$\frac{d\vec{M}_\perp}{dt} = -\vec{M}_\perp \wedge \vec{B}_0. \quad (8)$$

The last equation describes the precession of \vec{M}_\perp around \vec{B}_0 . The angular frequency of this process is called Larmor frequency

$$\omega_L = \gamma B_0. \quad (9)$$

Generating a transverse or anti-parallel magnetization to \vec{B}_0 can be achieved by applying a high frequency pulse to the ground state magnetization \vec{M}_0 . Let the static magnetic field \vec{B}_0 be pointing in z-direction. Now consider a coil oriented along the x-axis, if a sinusoidal voltage with frequency ω_H is applied it would generate a solenoidal magnetic field $\vec{B}_1(t)$ polarized along the x-direction. Under this conditions the torque acquires a second term:

$$\vec{\tau} = \vec{M}_0 \wedge (\vec{B}_0 + \vec{B}_1(t)) \quad (10)$$

which induces a second precession around the x-axis. For pulse durations which are short relative to the relaxation times, the magnitude of the magnetization is approximately constant. In this case the vector \vec{M} moves in

a sphere with radius $|\vec{M}_0|$. The vector coordinates in the sphere are then given by a azimuthal angle φ which is defined by the precession induced by \vec{B}_0 , and by a polar angle θ which arise from interaction with the solenoid field. Both angles are a functions of time:

$$\varphi = \omega_L t \quad (11)$$

$$\theta = \gamma B_1 t. \quad (12)$$

A pulse which induces a polar angle $\theta = 90^\circ$ is called a 90° pulse. For such a pulse the ground state magnetization is transformed into transverse magnetization \vec{M}_\perp , with $\vec{M}_\parallel = 0$. Analogously a pulse which generates $\theta = 180^\circ$ is called a 180° pulse. Such a pulse transforms the ground state into an anti-parallel magnetization $\vec{M}_\perp = -\vec{M}_0$, with $\vec{M}_\parallel = 0$.

1.2 NMR signal

1.2.1 Signal generation

In the set up used for this experiment the high frequency pulses of fixed frequency $\omega_{HF} = 19.8$ MHz were generated in a electronic unit of the minispec p20. Said pulses induce a precession around de z-axis, as explain in the previous section, that change the magnetic flux through the coil in time resulting in a induction signal modulated by the Larmor frequency ω_L , which could be set up by hand. This signal is fed back to the p20 electronic unit, where both the high frequency and the induction signals are mixed into an output given by the multiplication of both inputs, i.g the sum of two cosine functions. One of the terms depends on the working frequency, which is given by the difference between ω_L and ω_{HF} , and is on the order of few hundred Hertz while the second one in the range of 40MHz. The use of a low frequency bandpass filter allows to get rid of the former signal.

1.3 Relaxation Time

For the description of the precession of \vec{M}_\perp it is possible to define a rotating reference system (x', y', z) . This reference frame is defined by 2 condition: 1) the (x', y') -plane rotate in the static (x, y) -plane and 2) \vec{M}_\perp points in the x' -direction.

In this reference system the transverse and longitudinal magnetizations are constant if no relaxations processes are present, otherwise both components are time dependent. Their time evolution is described by the Bloch

equations:

$$\frac{d\vec{M}_{\perp}^{rot}}{dt} = -\frac{\vec{M}_{\perp}^{rot}}{T_2} \quad (13)$$

$$\frac{d\vec{M}_{\parallel}^{rot}}{dt} = -\frac{\vec{M}_{\parallel}^{rot} - \vec{M}^{rot}}{T_1}. \quad (14)$$

In the later equations it is assumed that the time evolution is dominated by a restoring force which is linearly proportional to the deflection from equilibrium. The constant T_2 in 13 is the so called spin-spin relaxation time whereas T_1 in 14 is called the spin-lattice relaxation time. The ground state is equal in the rotating and the static reference frame.

In the presence of relaxation processed the time evolution of the magnetization in laboratory \vec{M} and in rotating system \vec{M}^{rot} are related by

$$\frac{d\vec{M}}{dt} = \frac{\partial \vec{M}^{rot}}{\partial t} + \vec{\omega} \wedge \vec{M} \quad (15)$$

given the following Bloch equations in laboratory system:

$$\frac{d\vec{M}_{\perp}(t)}{dt} = \gamma(\vec{B} \wedge \vec{M})_{\perp} - \frac{\vec{M}_{\perp}(T)}{T_2}, \quad (16)$$

$$\frac{d\vec{M}_{\parallel}(t)}{dt} = \gamma(\vec{B} \wedge \vec{M})_{\parallel} - \frac{\vec{M}_{\parallel}(t) - \vec{M}_0}{T_1}. \quad (17)$$

1.4 Chemical shift

1.5 Imagin with NMR

2 Measurements

3 Analysis

4 Critical Discussion